

OXFORD

patient
handling



OXFORD

Voyager

Portable Overhead Lifter
and Easytrack™

*A lightweight ceiling hoist with
a unique portable tracking system.*

The options and accessories shown in the diagrams are available at additional cost(s).

 **SUNRISE
MEDICAL**



OXFORD

Voyager

Portable Overhead Lifter and Easytrack™

The Voyager Portable and Easytrack™ combine to form a unique portable overhead lifting solution for the home, care home and hospital. This combination of products have taken patient handling to the next dimension by providing the first truly portable overhead lifting solution. Whether it is for homecare, an emergency or temporary solution, or just that a fixed ceiling lifter is not suitable.... the Voyager Portable and Easytrack™ with a 200kg (31st) lifting capacity provides the perfect answer.

OXFORD

Voyager

Easytrack™

A unique overhead portable tracking system for the Voyager Portable hoist that requires no fixings, no alterations and can be installed within 5 minutes by a single person. The Easytrack™ is the perfect solution for homecare, temporary or emergency installations or even when the fixed Voyager Ceiling Track is not suitable.

Multi-Post Options – The 2, 3 or 4 post options allows the Easytrack™ to solve patient handling problems in nearly every room layout. The 2 post is ideal for a simple lift in or out of a bed, the 3 post for access in two directions and the 4 post for complete coverage.

Easy Set-up – The lightweight components allow any combination of the Easytrack™ to be set-up by a single person in under 5 minutes.

No Installation – The unique floor-to-ceiling design prevents any fixings, alterations or electrical installations allowing the system to be used permanently, semi-permanently or even temporarily. The Easytrack™ is ideal where the fixed Voyager Ceiling Track cannot be installed.

Portable Kit – The lightweight aluminium Easytrack™ comes in a portable kit form that can be easily transported along with the Voyager Portable hoist. A useful carry bag is available for both the Easytrack™ and Voyager Portable.

Footprint – The small feet of the Easytrack™ allow use even in the smallest bedroom or bathroom where a conventional floor lifter could not be manoeuvred.

Easy-Glide Rail – The unique design of the telescopic rail ensures that nearly any patient can be transferred without strain or risk of back injury.

Extending Rail – Unlike regular gantry tracks the Easytrack™ can be extended to suite different locations as and when they are required. The track can be extended when in position to allow adjustment and can reduce to ease storage or transport.

Telescopic Post – Allows the Easytrack™ to be used with a variety of ceiling heights.

Rotating Rail Connection – This allows the rail or post to be positioned at any angle, no matter how close to a wall or obstruction.

Patented Safety Systems – A multitude of unique safety systems have been designed into the Easytrack™ to ensure peace of mind. These systems include preventing the post losing grip if the floor or ceiling flexes, preventing the post being lowered if the rail is still connected or the rail disconnecting from the post.



S
e
t
u
p
i
n
4
e
a
s
y
s
t
e
p
s



Extend posts and rotate to lock into place



Fix one end of rail onto post.



Extend rail to required length and fix rail onto 2nd post.



Fix end cap to lock rail into place.

Easytrack™ Options

There are many variations of the Easytrack based around 2, 3 or 4 post models



The 2 post Easytrack™ with wall bracket



The 2 post Easytrack™ with bath bracket



The 2 post Easytrack™ without posts

The standard 2 post Easytrack™ is suited for transfer over a bed or bath in a peninsular setting.



The 3 post Easytrack™ with wall bracket



The 3 post Easytrack™ with bath bracket



The 3 post Easytrack™ without posts.

The 3 post Easytrack™ is suited for transfer over a bed or bath and down a room. Negates the need for a traditional curved track.



The 4 post Easytrack™ operates like a traditional xy system covering all 4 corners of a room. The 4 post is not telereafic.

An xy system without posts and with 4 wall brackets

OXFORD

Voyager

The lightweight
Voyager Portable
weighs only 5.5 kg (12 lbs)
but will lift a substantial
200 kg (440 lbs/31 st)



- 1 *Simple Trolley* allows the Voyager portable to be used with either the Easytrack™ or the standard Voyager Ceiling Track.
- 2 *Unique Quick Release Strap* allows instant repositioning of the hoist, saving time and battery power. Prevents ever having to lift the hoist onto the rail.
- 3 *Ergonomic Carry Handle* allows safe transportation of the hoist.
- 4 *Emergency Stop, Lower and strap twist prevention* helps ensure safety never becomes an issue.
- 5 *The Removable, Chargeable Battery* allows 24hr use and ensures the hoist is never out of action. An additional battery can be on charge when the other is in use.
- 6 *Secure Spreader Bar*, compatible with the full range of Oxford slings.
- 7 *Audible Low Battery Indicator* will help ensure the battery will not run-out at a critical moment.
- 8 *Rubber Feet* allow the hoist to stand upright safely and securely when stored.
- 9 *Soft Start and Stop* ensure the patient is transferred safely and comfortably.
- 10 *Dual Controls* coupled with a removable Handset allow total flexibility.
- 11 *Removable Handset* for added control.



Specifications

VOYAGER PORTABLE OVERHEAD HOIST		VOY-9800005.04-UK
Features		
Weight	5.5kg (12lbs) with battery	
Maximum User Weight	200kg (44lbs / 31st)	
Powered Vertically	4cm/sec	
Soft Start and Stop	✓	
Emergency Lower	Manual	
Emergency Stop	✓	
Removable Chargeable Battery	2 x 12v dc. Upto 30 transfers per charge	
Low Battery Indicator	Audible	
Freestanding Battery Charger	Holds one battery	
Quick Release Tape Mechanism	0 – 2.26m (88in). Manual release.	
Compatible with full Oxford Sling Range	✓	
Accessories		
	Code	
Additional Battery	VOY-A8500	
Additional Battery Charger	VOY-A8200	
Arm Extension	VOY-A8300	
Adapter Trolley for Voyager fixed Ceiling Track	VOY-A8900	

EASYTRACK™	2 Post System	3 Post System	4 Post System
Max User Weight	200 kg (440lbs/31st)	200 kg (440lbs/31st)	200 kg (440lbs/31st)
Ceiling Height Range	2130-2745mm (84-108")	2130-2745mm (84-108")	2130-2745mm (84-108")
Rail Length Range	1900-3160mm (74-124.5")	1900-3160mm (74-124.5")	3000mm (118")
Total Weight	20kg (44lbs)	42.5kg (94lbs)	57kg (125lbs)
Rail Weight (each)	12kg (26lbs)	12kg (26.4lbs)	9kg (20lbs)
Post Weight (each)	4kg (9lbs)	4kg (9lbs)	4kg (9lbs)
Standard Kit	VOY-9220010	VOY-9330010	VOY-9440030
Wall Bracket Option	VOY-9211010	VOY-9321010	-
Bath Bracket Option	VOY-9220110	VOY-9330110	-
Multi-Wall Bracket Option	VOY-9202010	VOY-9303010	VOY-9404030

Distributor:



Sunrise Medical Limited
High Street Wollaston
West Midlands DY8 4PS
Phone 01384 446688
Fax 01384 446699
www.sunrisemedical.co.uk

